

Accreditation Case Review: Mammography and Stereotactic Biopsy

Brett T. Parkinson MD

Breast Imaging Director Breast Care Services
Intermountain Healthcare

ACR Mammography Accreditation
Committee Chair

The Janice Beesley Hartvigsen Breast Care Center at Intermountain Medical Center



Clinical Image Evaluation

Scoring Procedures for Screen Film and
Full Field Digital Mammography

Mammography Accreditation Program

Accreditation Objectives

- Improve quality of mammography
- Reason for clinical image review
 - Quality of patient imaging is key
 - Provide feedback for improvement
- Uniform assessment of performance


Course Objectives

- Reasons for failed clinical images
 - Categories
 - Deficiencies
 - Problem solving and correcting deficiencies
 - Improving overall quality of images




Course Objectives

ACR Mammography Accreditation Program Pass Rates



Course Objectives


MAP Reasons for Unit Deficiency - 1ST ATTEMPT 2014



Course Objectives

Imaging Category	Failure Rate (%)
Positioning	1250 (20)
Exposure	944 (15)
Compression	887 (14)
Sharpness	806 (13)
Contrast	785 (13)
Artifacts	703 (11)
Labeling	465 (8)
Noise	288 (5)
Total	6128 (100)

Breast Density Past: The Old Descriptors



Fatty

Dense

Type 1	Type 2	Type 3	Type 4
Almost entirely fat (<10% fibroglandular)	Scattered fibroglandular densities (10%-50%)	Heterogeneously dense (51%-75% fibroglandular)	Extremely dense (>75% fibroglandular)


Breast Density Present: BI-RADS 2013

Breast Composition Categories


- The breasts are almost entirely fatty
- There are scattered areas of fibroglandular density
- The breasts are heterogeneously dense, which may obscure small masses
- The breasts are extremely dense, which lowers the sensitivity of mammography

(ACR BI-RADS ATLAS, 2013)


Fatty




Scattered



Heterogeneously Dense



Extremely Dense



Breast Density Future

We await publication of robust volume-based breast density data, using validated percentage cut points (not necessarily quartiles) that are readily and reproducibly determined at imaging, before again indicating percentage ranges for BI-RADS density categories.

(ACR BI-RADS ATLAS, p 126, 2013)

Overall Clinical Evaluation Form

ACR
American College of Radiology

Radiologist Clinical Image Evaluation
Facility MAP No. _____ Unit No. _____
IPA No. _____

Mammography Accreditation Program
PRIVILEGED and CONFIDENTIAL PEER REVIEW
Code of Virginia § 89.1-17

Reviewer No. _____ Date _____
Reviewer Name _____ Breast Type: FAT DENSE
 CLINICAL (Validation Film Check)

Justify any score of 1, 2, or 3. A score of 1 should only be used for the most serious problems and may require a review of additional images.

A. POSITIONING

	RCC	LCC	RMLO	LMLO
<input type="checkbox"/> 1 Poor visualization of posterior tissues				
<input type="checkbox"/> 2 Sagging breast				
<input type="checkbox"/> 3 Inadequate amount of pectoral muscle shown on image				
<input type="checkbox"/> 4 Inadequate inframammary fold (IMF)				
<input type="checkbox"/> 5 Poor visualization of posterior tissues				
<input type="checkbox"/> 6 Excessive exaggeration				
<input type="checkbox"/> 7 Portion of breast cut off				
<input type="checkbox"/> 8 Skin folds				
<input type="checkbox"/> 9 Other body parts projected over breast				
<input type="checkbox"/> 10 Breast positioned too high on image				
<input type="checkbox"/> 11 Posterior nipple line (PNL) on CC not within 1 cm of MLO PNL				

1 Technologist technique 5 Uncertain
 2 Inappropriate mammographic projections
 3 Wrong size image receptor
 4 Other: _____

Additional comments on positioning: _____

B. COMPRESSION

	RCC	LCC	RMLO	LMLO
<input type="checkbox"/> 1 Poor separation of parenchymal densities				
<input type="checkbox"/> 2 Non-uniform exposure levels				
<input type="checkbox"/> 3 Patient motion				
<input type="checkbox"/> 4 Other				

1 Under compression by technologist
 2 Unsuitable compression device
 3 Technologist positioning of compression device
 4 Other: _____
 5 Uncertain

C. EXPOSURE LEVEL

	RCC	LCC	RMLO	LMLO
<input type="checkbox"/> 1 Generalized underexposure				
<input type="checkbox"/> 2 Generalized overexposure				
<input type="checkbox"/> 3 Inadequate penetration of dense areas				
<input type="checkbox"/> 4 Excessive penetration of lucent areas				
<input type="checkbox"/> 5 Other				

1 Film development 7 Uncertain
 2 Under compression with phototiming
 3 Radiologist preference
 4 Phototimer variability
 5 Other: _____

D. CONTRAST

	RCC	LCC	RMLO	LMLO
<input type="checkbox"/> 1 Inadequate contrast				
<input type="checkbox"/> 2 Excessive contrast				
<input type="checkbox"/> 3 Other				


1 Film development 2 Improper kVp
 3 Excessive scatter 4 Uncertain
 4 Underexposure
 5 Digital: window width too wide
 6 Digital: window width too narrow
 7 Other: _____

Deficiencies	BreastView	Possible Causes	Score				
E. SHARPNESS			1 2 3 4 5				
<input type="checkbox"/> 1 Poor delineation of linear structures	RCC LCC RMLO LMLO	<input type="checkbox"/> Patient motion	<input type="checkbox"/> 5 Uncertain				
<input type="checkbox"/> 2 Poor delineation of feature margins		<input type="checkbox"/> Poor screen contact					
<input type="checkbox"/> 3 Poor delineation of microcalcifications		<input type="checkbox"/> Film-screen separation					
		<input type="checkbox"/> Other: _____					
F. NOISE			1 2 3 4 5				
<input type="checkbox"/> 1 Visually striking mottle pattern	RCC LCC RMLO LMLO	<input type="checkbox"/> Film development	<input type="checkbox"/> 3 Improper kVp				
<input type="checkbox"/> 2 Noise limited visualization of detail		<input type="checkbox"/> Recording system speed	<input type="checkbox"/> 5 Uncertain				
		<input type="checkbox"/> Digital: inadequate SNR					
		<input type="checkbox"/> Digital: window width too narrow					
		<input type="checkbox"/> Other: _____					
G. ARTIFACTS			1 2 3 4 5				
<input type="checkbox"/> 1 Punctate or lint	LCC RCC LMLO RMLO	<input type="checkbox"/> Poor screen maintenance	<input type="checkbox"/> 5 Uncertain				
<input type="checkbox"/> 2 Scratches or pickoff		<input type="checkbox"/> Development related					
<input type="checkbox"/> 3 Roller marks		<input type="checkbox"/> Unsuitable grid or Bucky					
<input type="checkbox"/> 4 Grid related artifacts		<input type="checkbox"/> Film exposed to light					
<input type="checkbox"/> 5 Hair, deodorant, etc.		<input type="checkbox"/> Lack of patient preparation					
<input type="checkbox"/> 6 Film handling		<input type="checkbox"/> Poor cassette closure					
<input type="checkbox"/> 7 Film fogging		<input type="checkbox"/> Damaged cassette					
<input type="checkbox"/> 8 Poor screen-film alignment		<input type="checkbox"/> Digital: detector calibration (eg, uniformity calibration)					
<input type="checkbox"/> 9 Digital: image receptor artifact		<input type="checkbox"/> Digital: foreign objects calibrated into calibration file					
<input type="checkbox"/> 10 Digital: laser printer artifact		<input type="checkbox"/> Digital: laser printer needs service					
<input type="checkbox"/> 11 Digital: laser printer scanning lines		<input type="checkbox"/> Other: _____					
<input type="checkbox"/> 12 Other							
H. EXAM ID (Only fail this category if the patient's name, ID or laterality is missing)			1 2 3 4 5				
<input type="checkbox"/> 1 Patient name and additional patient identifier	RCC LCC RMLO LMLO	<input type="checkbox"/> Technologist error	<input type="checkbox"/> 5 Uncertain				
<input type="checkbox"/> 2 Facility name and location (city, state and zip)		<input type="checkbox"/> Missing or non-standard labeling method					
<input type="checkbox"/> 3 Date of examination		<input type="checkbox"/> Improper positioning of label					
<input type="checkbox"/> 4 View and laterality		<input type="checkbox"/> Other: _____					
<input type="checkbox"/> 5 Unit identification (if more than one)							
<input type="checkbox"/> 6 Technologist identification							
<input type="checkbox"/> 7 Cassette/screen identification							
OVERALL CLINICAL IMAGE EVALUATION RECOMMENDATION			<input type="checkbox"/> 2 Pass <input type="checkbox"/> 1 Fail				
Circle categories leading to a failing grade. Only fail for artifacts or exam ID problems if they are severe.							
A	B	C	D	E	F	G	H
Positioning	Compression	Exposure	Contrast	Sharpness	Noise	Artifacts	Exam ID
ADDITIONAL RECOMMENDATIONS:							
<input type="checkbox"/> 1 Recommend reviewing 1999 ACR QC Manual for MQSA film labeling requirements							
<input type="checkbox"/> 2 Recommend sending technologist to a hands-on positioning course							
<input type="checkbox"/> 3 Excessive collimation: recommend coning to size of film to reduce viewbox glare							
<input type="checkbox"/> 4 Check optimal film development conditions with film manufacturer							
<input type="checkbox"/> 5 Physician should review Clinical Image Evaluation section of 1999 ACR QC Manual							
ADDITIONAL COMMENTS: (If you find an abnormality, mark it with a grease pencil and comment below)							

Positioning: MLO Most Common Deficiencies

Deficiency	Frequency (%)
Inadequate pectoral muscle	35
“Sagging” of the breast	22
Poor visualization posterior tissue	22
Skin folds overlying breast tissue	10
Breast positioned too high on image receptor	6
Portion of breast cut off	5


Positioning: MLO Most Common Deficiencies



FG = fibroglandular tissue; M = Pectoralis muscle; N = Nipple;
IF = inframammary fold


Adequate Pectoral Muscle

RMLO




- Optimal pectoral muscle on a digital image
- Length
 - To nipple line or below
- Shape
 - Convex anterior border
 - Wider superiorly
 - Gradually narrows inferiorly
- Score = 5

Inadequate Pectoral Muscle



- When compared with a straight line, at a minimum the muscle should align with the line or go beyond
- Preferably, it will bulge beyond the edge
- Note difference in width of superior vs inferior aspect of muscle
- Score = 4

Sagging Breast: Poor Delineation of Structures



Inadequate muscle


Overlapping breast tissues

N

skin fold

Score = 2


To prevent sagging: Out and Nipple-Up Maneuver



To prevent sagging: Out and Nipple-Up Maneuver

Sagging breast


Repeated exam




skin fold

Score = 5

Sagging Breasts with Inframammary Skin Folds




Open IMF



Open Inframammary Fold

- Skin overlying abdomen in front of edge of receptor
- Free of skin folds



Positioning Criteria: CC

- Inclusion of medial tissue
 - Nipple centered on film
 - No excessive exaggeration
 - Exaggeration may cause loss of medial or lateral posterior tissue
 - Inclusion of posterior tissue
 - Free of skin folds
-


Positioning: CC Most Common Deficiencies



Deficiency	Frequency (%)
Poor visualization of posterior tissue on CC	32
Posterior nipple line on CC not within 1 cm of MLO	30
Excessive lateral/medial exaggeration on CC	20
Skin folds overlying breast tissue	18

Visualization of Posterior Breast on CC



- PNL should measure within 1 cm of MLO
 - Depends on a well positioned MLO
- Pectoral muscle can be seen in 30% of women
- Requires proper positioning of CC:
 - Elevate inferior breast (freely movable tissue)
 - Pull superior and inferior tissue onto receptor
 - Lean patient's head forward to the side of tube
- Pectoral muscle in central and medial breast
 - Improper positioning if central and lateral



Visualization of Posterior Breast on CC




Radiologist Clinical Image Evaluation
Facility MAP No. _____ Unit No. _____
IPA No. _____



Inclusion of Medial Tissue


- Skin reflection of cleavage visualized
- Retroglandular fat
 - Sternum flush with bucky
 - Contralateral breast draped
- Free of skin folds
 - IMF properly elevated
 - Check for gap
 - Must look from medial side



Nipple Centered on Film

- No excessive exaggeration
- Otherwise, nipple will point to the missing tissue
 - If points laterally RT should lift and pull out lateral tissue
 - If points medially RT should rotate to place sternum flush with bucky


Score = 4





Excessive Exaggeration

- Nipple will point to the laterality of missing tissue
 - If it points laterally, lift and pull out lateral tissue
 - If points medially, rotate to place sternum flush with bucky

Score = 3




Skin Folds Overlying Breast



Score = 2


Inframammary skin folds



Score = 3 or 2?

Depend on other images; Pec here = 2

Medial skin folds on CC



Compression


- Affects sharpness
 - Compression decreases unsharpness (blur) due to reduced time and motion
- Affects contrast
 - Good compression increases contrast due to less scattered radiation
- Affects exposure level
 - Compression increases exposure level due to reduced breast thickness

Compression: Most Common Deficiencies

Deficiency	Frequency (%)
Poor separation of fibroglandular tissues	65
Non-uniform exposure of fibroglandular tissues	19
Patient motion	16

Inadequate Compression


- Poor separation of fibroglandular tissues
- Unequal exposure of fibroglandular tissues
- Motion unsharpness




Assessing Sharpness on the MLO

- Anterior edge of pectoral muscle (blurred?)
- Subareolar area
- Inferior breast
 - Edges of vessels, calcifications, and Cooper's ligaments


Score = 2




Compression



Year 1
Score = 2




Year 2
Score = 2-3



Year 3
Score = 4

Assessing Sharpness on CC

- Subareolar area
- Medial breast
 - Edges of vessels, calcification, ligaments




Assessing Sharpness on CC

- Minimum average optical density (OD) in fibroglandular tissue should be 1.0 when measured with a spot densitometer
- Fibroglandular tissue should be a shade of gray
- Pectoral muscle may have OD under 1.0; however, must be able to see underlying tissue
- Uniformly washed out look indicates under- exposure on screen film images

Exposure Level: Most Common Deficiencies

Deficiency	Frequency (%)	
Generalized underexposure	56	} (80)
Inadequate exposure of dense tissues	24	
Generalized overexposure	18	} (20)
Overexposure of radiolucent tissues	2	

Exposure Level: Most Common Deficiencies




Contrast

Differences in optical densities between different tissues. Fat should have high optical density (black) and fibro-glandular tissue should have much lower optical density (white).


Contrast: Most Common Deficiencies

Deficiency	Frequency (%)
Inadequate contrast	90
Excessive contrast	10


Contrast



Inadequate




Good



Excessive

Good Contrast Film-Screen and Digital: Same Patient



Sharpness


- Unsharpness (often referred to as “blur” by medical physicists)
 - Blurring of margins of vessels and Cooper’s ligaments
 - Unclear margins of masses
 - Unsharpness of calcifications
-

Assessing Sharpness on the MLO



- Inferior third of the breast
- Subareolar area
- Anterior edge of pectoral muscle
- Edges of vessels, calcification, and Cooper’s ligaments

Inadequate Compression

- Motion unsharpness seen as poor delineation of linear structures
- This is usually most evident in inferior and subareolar aspect of MLO
- Also look at the edge of the pectoral muscle



Inadequate Compression




Motion unsharpness

Score = 3

Assessing Sharpness on the CC

- Medial aspect of the breast
 - Subareolar area
 - Edges of vessels, calcifications, Cooper's ligament
-

Assessing Sharpness on the CC




- Unsharpness related to inadequate compression and motion more likely to occur on MLO
- If question about sharpness, compare CC to MLO

Noise

- Mottled pattern of areas of relatively equal tissue density
 - Noise limits visualization of details in the image
-


Noise

Mottled pattern of areas of relatively equal tissue density



Noise

If noise suspected, it helps to use a mag lens



Noise

Mottled pattern makes it difficult to see fine details in the images



Scoring Categories

Any score 3 or below - reviewer must indicate:

- Deficiency(s)
- Possible cause(s)
- And provide comments, if possible

Justify any score of 1, 2 or 3. A score of 1 should **only** be used for the most serious problems and may require a review of additional images.

Deficiencies	Breast/ View				Possible Causes	Score
	RCC	LCC	RMLO	LMLO		
A. POSITIONING						WORST (circle ONLY one) BEST 1 2 3 4 5
<input type="checkbox"/> ¹ Poor visualization of posterior tissues					<input type="checkbox"/> ¹ Technologist technique	<input type="checkbox"/> ⁵ Uncertain
<input type="checkbox"/> ² Sagging breast					<input type="checkbox"/> ² Inappropriate mammographic projections	
<input type="checkbox"/> ³ Inadequate amount of pectoral muscle shown on image					<input type="checkbox"/> ³ Wrong size image receptor	
<input type="checkbox"/> ¹³ Inadequate inframammary fold (IMF)					<input type="checkbox"/> ⁴ Other: _____	
<input type="checkbox"/> ⁴ Poor visualization of posterior tissues					Additional comments on positioning:	
<input type="checkbox"/> ⁵ Excessive exaggeration					_____	
<input type="checkbox"/> ⁶ Portion of breast cut off					_____	
<input type="checkbox"/> ⁷ Skin folds					_____	
<input type="checkbox"/> ⁸ Other body parts projected over breast					_____	
<input type="checkbox"/> ¹² Breast positioned too high on image receptor					_____	
<input type="checkbox"/> ¹¹ Posterior nipple line (PNL) on CC <i>not</i> within 1 cm of MLO PNL					_____	

Failing Score

- Score of 1 or 2 in any category
- 3 (or more) scores of 3 in categories (A – G)
 - Does *not* include category H (Exam ID)
- Do *NOT* fail for minor artifacts (G)
 - Artifacts should be prominent to earn a score of 3
 - Artifacts should be severe to earn a score of 2 or 1
- Do fail H (exam ID) even with a score of 2 (the patient's name, patient ID or R/L are missing)
 - Do *NOT* fail for other minor exam ID deficiencies

Case 1

